



Anti-ZBTB2 (aa 50-100) polyclonal antibody (CABT-BL3831)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within amino acids 50-100 of Human ZBTB2 according to NP_065912.1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Cellular Localization	Nucleus.
Format	Liquid
Buffer	99% Tris citrate/phosphateNote: pH: 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	ZBTB2 (zinc finger and BTB domain containing 2) is a protein-coding gene. Diseases associated with ZBTB2 include chronic myeloid leukemia, and myeloid leukemia. GO annotations related to this gene include DNA binding and zinc ion binding. An important paralog
---------------------	---

of this gene is ZBTB25.

GENE INFORMATION

Entrez Gene ID	57621
----------------	-----------------------

Protein Refseq	NP_065912
----------------	---------------------------

UniProt ID	Q8N680
------------	------------------------
