



Anti-ZDHHC4 polyclonal antibody (CABT-BL3836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human ZDHHC4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Probable palmitoyltransferase ZDHHC4 is an enzyme that in humans is encoded by the ZDHHC4 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID [55146](#)

Protein Refseq [NP_001127859](#)

UniProt ID [A0A024QZS8](#)
