



Anti-ZFP37 (aa 252-301) polyclonal antibody (CABT-BL3838)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 252-301 (LCCHSSSHIK QDKIQTGEKH EKSPSLSSST KHEKPQACVK PYECNQCGKV) of Human ZFP37.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	ZFP37 may be involved in transcriptional regulation.
---------------------	--

GENE INFORMATION

Entrez Gene ID [7539](#)

Protein Refseq [NP_003399](#)

UniProt ID [Q9Y6Q3](#)
