



# Anti-ZFP472 (aa 20-69) polyclonal antibody (CABT-BL3840)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within N terminal amino acids 20-69 (TLGEWTMLDS SQKKLYRDVM KETFLNLISI GKTEEDTEE EYQNPKNLR) of Mouse Zfp472.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	97% PBS, 2% Sucrose
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">224691</a>
<b>Protein Refseq</b>	<a href="#">NP_694703</a>

UniProt ID

[B0V2W5](#)

---