



Anti-ZFP759 (aa 30-79) polyclonal antibody (CABT-BL3841)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 30-79 (LGPPQWKLYR DVMLENNYL VFLGLASSKP YLVRFLAQIQ EPLDVKSQVA) of Mouse Zfp759 (EDL10320.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	268670
Protein Refseq	NP_765980

UniProt ID

[Q7M6X3](#)
