



Anti-ZFY1 (aa 326-375) polyclonal antibody (CABT-BL3842)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: VSEIKAAFLP IAWTAAYDNN SDEIEDQNVT ASALLNQDES GGLDRVPKQK, corresponding to internal sequence amino acids 326-375 of Mouse ZFY
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	22767
Protein Refseq	NP_033596

UniProt ID

[P10925](#)
