



Anti-ZK1251.1 (aa 1-5/81-96) polyclonal antibody (CABT-BL3847)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptides SGRGK and RHLQLAIRNDEELNKL, corresponding to amino acids 1-5 and 81-96 of Human Histone H2A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, <i>S. cerevisiae</i> , <i>C. elegans</i> , Fruit fly, <i>S. pombe</i>
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ICC, IHC-P, WB, ChIP
Cellular Localization	Nucleus. Chromosome.
Format	Liquid
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [191549](#)

Protein Refseq [NP_501718](#)

UniProt ID [Q23429](#)
