



Anti-ZKSCAN14 (aa 301-350) polyclonal antibody (CABT-BL3848)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 301-350 (SVKPHSSVDS AVGLLETQRQ FQEDKPYKCD SCEKGFRQRS DLFKHQRIHT) of Mouse ZNF394
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [67235](#)

Protein Refseq [NP_075811](#)

UniProt ID [Q9Z1D9](#)
