



Anti-ZKSCAN2 (aa 917-967) polyclonal antibody (CABT-BL3849)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (Human) corresponding to a region between residues 917 and 967 of ZKSCAN2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	99% Tris citrate/phosphateNote: pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	342357
Protein Refseq	NP_001012999

UniProt ID

[Q63HK3](#)
