



Anti-ZMYND15 polyclonal antibody (CABT-BL3850)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human ZMYND15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	ZMYND15 (zinc finger, MYND-type containing 15) is a protein-coding gene. GO annotations related to this gene include histone deacetylase binding and zinc ion binding.
--------------	--

GENE INFORMATION

Entrez Gene ID [84225](#)

Protein Refseq [NP_115641.1](#)

UniProt ID [Q9H091](#)
