



Anti-ZMYND17 (aa 326-375) polyclonal antibody (CABT-BL3851)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 326-375 (SAHRGLYHDF WEEQVETGQT HHPDLVAAFH PGFHSSPDLM EAWLPTLLLL) of Human ZMYND17.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	118490
Protein Refseq	NP_001019764

