



Anti-ZNF132 (aa 215-264) polyclonal antibody (CABT-BL3856)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal sequence amino acids 215-264 (GQKPYSNLGQ LPEVCTTQKL FECSNCGKAF LKSSTLPNHL RTHSEEIPFT) of Human ZNF132, NP_003424.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7691](#)

Protein Refseq [NP_003424](#)

UniProt ID [B3KQ54](#)
