



Anti-ZNF135 (aa 29-78) polyclonal antibody (CABT-BL3858)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within a region corresponding to N terminal amino acids 29-78 (ELWAVESRLP QGVYPDLETR PKVKLSVLKQ GISEEISNSV ILVERFLWDG) of Human ZNF135.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7694](#)

Protein Refseq [NP_009065](#)

UniProt ID [Q8N9M3](#)
