



Anti-ZNF138 (aa 215-262) polyclonal antibody (CABT-BL3860)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|--|
| Immunogen | Synthetic peptide corresponding to a region within C terminal amino acids 215-262 of Human zinc finger protein 138, IHTTEKPYKC EQCGKVKQSQ PTLTKHQIY TGEEPYKCEE CGKAFNLS. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB |
| Cellular Localization | Nuclear |
| Format | Liquid |
| Buffer | 2% Sucrose, PBS |
| Preservative | None |
| Storage | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles. |

GENE INFORMATION

Entrez Gene ID [7697](#)

Protein Refseq [NP_001153655](#)

UniProt ID [A2RRP7](#)
