



Anti-ZNF138 (aa 1-172) polyclonal antibody (CABT-BL3861)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to N terminal amino acids 1-172 of Human zinc finger protein 138 (BC131762).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	98% PBS, 1% BSA pH: 7.20
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	7697
Protein Refseq	NP_001153655

UniProt ID

[A2RRP7](#)
