



Anti-ZNF141 (aa 145-194) polyclonal antibody (CABT-BL3863)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: LQCKASVKVV SKFSNSNKRK TRHTGEKHFK ECGKSFQKFS HLTQHKVIHA, corresponding to N terminal amino acids 145-194 of Human ZNF141
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7700](#)

Protein Refseq [NP_003432](#)

UniProt ID [Q15928](#)
