



Anti-ZNF141 (aa 424-473) polyclonal antibody (CABT-BL3864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within internal sequence amino acids 424-473 (CKECDKAFKQ FSLLSQHKKI HTVDKPYKCK DCDKAFKRFS HLNKHKKIHT) of Human ZNF141
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Reconstitution	Add 50 µl of distilled water. Final antibody concentration is 1 mg/ml.
Cellular Localization	Nuclear
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7700](#)

Protein Refseq [NP_003432](#)

UniProt ID [Q15928](#)
