



Anti-ZNF157 (aa 2-51) polyclonal antibody (CABT-BL3870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 2-51 (PANGTSPQRF PALIPGEPGR SFEGSVSFED VAVDFTRQEW HRLDPAQRTM) of human ZNF157.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	May be involved in transcriptional regulation.
---------------------	--

GENE INFORMATION

Entrez Gene ID [7712](#)

Protein Refseq [NP_003437](#)

UniProt ID [P51786](#)
