



# Anti-ZNF17 polyclonal antibody (CABT-BL3872)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to a region of Human ZNF17
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID	<a href="#">7565</a>
Protein Refseq	<a href="#">NP_008890</a>

UniProt ID

[P17021](#)

---