



Anti-ZNF192 (aa 2-51) polyclonal antibody (CABT-BL3875)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within residues: AEESRKPSAP SPPDQTPEED LVIVKVEEDH GWDQESSLHE SNPLGQEVFR, corresponding to N terminal amino acids 2-51 of Human ZNF192
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Reconstitution	Add 50 µl of distilled water for a final concentration of 1 mg/ml.
Cellular Localization	Nucleus.
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7745](#)

Protein Refseq [NP_006289](#)

UniProt ID [A0A024RCK7](#)
