



Anti-ZNF219 (aa 161-191) polyclonal antibody (CABT-BL3881)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to N terminal amino acids 161-191 of Human Zfp219
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Ammonium Sulphate Precipitation
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Nucleus.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	May be involved in transcriptional regulation.
---------------------	--

GENE INFORMATION

Entrez Gene ID [51222](#)

Protein Refseq [NP_001095924](#)

UniProt ID [Q9P2Y4](#)
