



Anti-ZNF221 (aa 108-157) polyclonal antibody (CABT-BL3882)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 108-157 (METVPEAGPH EEWSCQIWE QIASDLTRSQ NSIRNSSQFF KEGDVPCQIE) of Human ZNF221.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7638](#)

Protein Refseq [NP_037491](#)

UniProt ID [A0A087WT08](#)
