



# Anti-ZNF260 polyclonal antibody (CABT-BL3886)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nucleus
Format	Lyophilized Powder
Buffer	PBS with 0.02% Sodium Azide, 1%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20°C, Avoid freeze / thaw cycles

## BACKGROUND

Introduction	Transcription factor that acts as a cardiac regulator and an effector of alpha1-adrenergic signaling. Binds to PE response elements (PERE) present in the promoter of genes such as ANF/NPPA and acts as a direct transcriptional activator of NPPA. Also acts as a cofactor with GATA4, a key cardiac regulator.
--------------	---

# GENE INFORMATION

Entrez Gene ID [339324](#)

---

Protein Refseq [NP\\_001159510.1](#)

---

UniProt ID [Q3ZCT1](#)

---