



Anti-ZNF281 (aa 846-895) polyclonal antibody (CABT-BL3891)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 846-895 (PMTFITNSNG EVDHRVRTSV SDFSGYTMM SDVSEPCSTR VKTPTSQSYR) of human Zfp281.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Involved in transcriptional regulation. Represses the transcription of a number of genes including gastrin and ornithine decarboxylase. Binds to the G-rich box in the enhancer region of
---------------------	---

these genes.

GENE INFORMATION

Entrez Gene ID	23528
----------------	-----------------------

Protein Refseq	NP_036614
----------------	---------------------------

UniProt ID	Q9Y2X9
------------	------------------------
