



Anti-ZNF295 polyclonal antibody (CABT-BL3894)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to a region of Human ZN295.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Zinc finger protein 295 is a protein that in humans is encoded by the ZNF295 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [49854](#)

Protein Refseq [NP_001091872](#)

UniProt ID [Q9ULJ3](#)
