



Anti-ZNF3 polyclonal antibody (CABT-BL3895)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region of human ZNF3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Zinc finger protein 3 is a protein that in humans is encoded by the ZNF3 gene.

GENE INFORMATION

Entrez Gene ID [7551](#)

Protein Refseq [NP_060185](#)

UniProt ID [P17036](#)
