



Anti-ZNF318 polyclonal antibody (CABT-BL3898)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human ZNF318.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Zinc finger protein 318 is a protein that in humans is encoded by the ZNF318 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [24149](#)

Protein Refseq [NP_055160](#)

UniProt ID [Q5VUA4](#)
