



# Anti-ZNF395 (internal region) polyclonal antibody (CABT-BL3905)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide (Human) from the center region of Papillomavirus Regulatory Factor 1 / HDBP2, conjugated to KLH.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, IHC-P
<b>Cellular Localization</b>	Cytoplasmic and Nuclear. Note= May shuttle between nucleus and cytoplasm.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Plays a role in papillomavirus genes transcription.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [55893](#)

---

Protein Refseq [NP\\_061130](#)

---

UniProt ID [Q7L9C8](#)

---