



# Anti-ZNF409 polyclonal antibody (CABT-BL3906)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide (Human). More details of the immunogen are not available.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P
<b>Cellular Localization</b>	Nuclear
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [22830](#)