



Anti-ZNF431 (aa 2-51) polyclonal antibody (CABT-BL3908)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N-terminal amino acids 2 - 51 (DDLKYGVYPLKEASGCPGAERNLLVSYFEKETLTFRDV AIEFSLEEWEC) of Human ZNF431
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [170959](#)

Protein Refseq [NP_597730](#)

UniProt ID [A0A024R7Q8](#)
