



Anti-ZNF440 (aa 38-87) polyclonal antibody (CABT-BL3910)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: RNLTSLGKRW KDQNIYEYHQ NPRRNFRSLI EEKVNEIKDD SHCGETFTPV, corresponding to amino acids 38-87 of Human ZNF440
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Reconstitution	Add 100 µl of distilled water.
Cellular Localization	Nuclear
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [126070](#)

Protein Refseq [NP_689570](#)

UniProt ID [Q8IYI8](#)
