



Anti-ZNF454 (aa 108-157) polyclonal antibody (CABT-BL3918)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 108-157 (WTIKERFSSS SHWKASLLE WQCGGQEISL QRVVLTHPNT PSQEYDESGS) of Human ZNF454, NP_872400.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [285676](#)

Protein Refseq [NP_872400](#)

UniProt ID [Q8N9F8](#)
