



Anti-ZNF474 (aa 252-301) polyclonal antibody (CABT-BL3919)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 252-301 (DRLPVELHQP LPQKPQPLPN AQSSQAGPNQ AQLVFCPHCS RIFTSDRLLV) of Human FLJ32921.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	133923
Protein Refseq	NP_997200

UniProt ID

[Q6S9Z5](#)
