



Anti-ZNF491 (aa 2-51) polyclonal antibody (CABT-BL3921)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 2-51 (GERLFESAEG SQCGETFTQV PEDMLNKKTL PGVKSCESGT CGEIFMGYSS) of human ZNF491.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	126069
Protein Refseq	NP_689569

UniProt ID

[Q8N8L2](#)
