



Anti-ZNF527 (aa 528-577) polyclonal antibody (CABT-BL3927)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: GKAFHQILSL RLHQRIHAGE KPYKCNECGN NFSCVSALRR HQRIHNRETL, corresponding to C terminal amino acids 528-577 of Human ZNF527
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [84503](#)

Protein Refseq [NP_115829](#)

UniProt ID [Q8NB42](#)
