



Anti-ZNF541 (aa 360-409) polyclonal antibody (CABT-BL3928)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A 14-22 amino acid region within synthetic peptide: ASPSQVAMAS FSSAGPPADP SKSKLTIFSR IQGGNIYRLP HPVKEENVAG, corresponding to internal sequence amino acids 360-409 of Human ZNF541
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [84215](#)

Protein Refseq [NP_001094889](#)

UniProt ID [Q9H0D2](#)
