



Anti-ZNF579 (aa 512-562) polyclonal antibody (CABT-BL3936)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within amino acids 512-562 of Human ZNF579.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Nucleus.
Format	Liquid
Buffer	99% Tris citrate/phosphateNote: pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	163033
Protein Refseq	NP_689813

UniProt ID

[Q8NAF0](#)
