



Anti-ZNF579 (aa 503-552) polyclonal antibody (CABT-BL3937)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal sequence amino acids 503-552 (LANIKEEPPS PGTPPQSPPA PPVFLSASCF DSQDHSAFEM EEEEVDSKAH) of Human ZNF579.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Reconstitution	Add 50 µl of distilled water for a final antibody concentration of 1 mg/ml.
Cellular Localization	Nucleus.
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [163033](#)

Protein Refseq [NP_689813](#)

UniProt ID [Q8NAF0](#)
