



Anti-ZNF583 (aa 72-121) polyclonal antibody (CABT-BL3938)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|--|
| Immunogen | Synthetic peptide from a region within N terminal amino acids 72-121 (TRGPCPDWEY VFKNSEFSSK QETYEESKV VTVGARHLSY SLDYPSLRED) of Human ZNF583 |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB |
| Reconstitution | Add 50 µl of distilled water for a final antibody concentration of 1 mg/ml. |
| Cellular Localization | Nuclear |
| Format | Lyophilised |
| Buffer | 2% Sucrose, PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|---------------------|--|
| Introduction | May be involved in transcriptional regulation. |
|---------------------|--|

GENE INFORMATION

Entrez Gene ID [147949](#)

Protein Refseq [NP_001153332](#)

UniProt ID [Q96ND8](#)
