



Anti-ZNF600 (internal region) polyclonal antibody (CABT-BL3939)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthesized peptide derived from an internal region of Human ZNF600.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nucleus.
Format	Liquid
Buffer	PBS, 0.88% Sodium chloride, 50% Glycerol pH: 7.40 Note: PBS (without Mg2+ and Ca2+)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	162966
Protein Refseq	NP_940859

UniProt ID

[Q6ZNG1](#)
