



Anti-ZNF609 (aa 72-121) polyclonal antibody (CABT-BL3940)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: IKFVTPVPGP QGKEGKSKSK RSKSGKDTSK PTPGTSLFTP SEGAASKKEV, corresponding to amino acids 72-121 of Human ZNF609
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus (Probable).
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [23060](#)

Protein Refseq [NP_055857](#)

UniProt ID [O15014](#)
