



Anti-ZNF627 (aa 37-86) polyclonal antibody (CABT-BL3946)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: FRNLASVGKQ WEDQNIEDPF KIPRRNISHI PERLCESKEG GQGEETFSQI, corresponding to a region within N terminal amino acids 37-86 of Human ZNF627
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	199692
-----------------------	------------------------

Protein Refseq [NP_660338](#)

UniProt ID [Q7L945](#)
