



Anti-ZNF645 polyclonal antibody (CABT-BL3952)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human ZNF645.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene encodes a member of the zinc finger domain-containing protein family. This family member contains both a RING-type and a C2H2-type of zinc finger domain, and it may function as an E3 ubiquitin-protein ligase. Protein localization suggests a role in human sperm production and quality control.
--------------	--

GENE INFORMATION

Entrez Gene ID [158506](#)

Protein Refseq [NP_689790](#)

UniProt ID [Q8N7E2](#)
