



Anti-ZNF646 (aa 218-267) polyclonal antibody (CABT-BL3953)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 218-267 (VVNFTGGQEP TQSPPAEEER RYKCSQCGKT YKHAGSLTNH RQSHTLGIYP) of Human ZNF646.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [9726](#)

Protein Refseq [NP_055514](#)

UniProt ID [O15015](#)
