



Anti-ZNF660 (aa 2-51) polyclonal antibody (CABT-BL3956)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 2-51 (RRKTRNFKHK TVKDNKVLTE GSDQSEKDN SQCCDPATNE RVQAEKRQYV) of Human ZNF660, NP_775929.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [285349](#)

Protein Refseq [NP_775929](#)

UniProt ID [Q6AZW8](#)
