



Anti-ZNF700 (aa 649-698) polyclonal antibody (CABT-BL3961)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the C terminal amino acids 649-698 (CKECEKAFCK FSSFQIHERK HRGEKPYECK HCGNGFTSAK ILQIHARTHI) of Human ZNF700.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	90592
-----------------------	-----------------------

Protein Refseq [NP_653167](#)

UniProt ID [A0A087WVH9](#)
