



Anti-ZNF708 (aa 433-482) polyclonal antibody (CABT-BL3963)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within C terminal amino acids 433-482 (ILTKHKVIHTEDKPYKCEECKGKTFNYSSNFTNHKKIHTG EKPYKCEECKGK) of Human ZNF708, EAW84891
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [7562](#)

Protein Refseq [NP_067092](#)

UniProt ID [P17019](#)
