



Anti-ZNF720 (aa 34-62) polyclonal antibody (CABT-BL3965)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide designed between amino acids 34-62 from the central region of Human ZNF720.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	124411
Protein Refseq	NP_001124385
UniProt ID	Q7Z2F6