



Anti-ZNF799 (aa 108-157) polyclonal antibody (CABT-BL3974)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: SGEVIMGHSS LNCYIRVGAG HKPYEYHECG EKPDTHKQRG KAFSYHNSLQ, corresponding to internal sequence amino acids 108-157 of Human ZNF799
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	90576
Protein Refseq	NP_001074290

UniProt ID

[Q96GE5](#)
