



Anti-ZNF823 (internal region) polyclonal antibody (CABT-BL3976)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A KLH conjugated synthetic peptide selected from the center region of Human HSZFP36.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Ammonium Sulphate Precipitation
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	55552
Protein Refseq	NP_001073962

UniProt ID

[P16415](#)
